

# Tracing metals and comparing BC sources in atmospheric aerosols samples collected at UALR

---

Angie Marchany-Rivera  
Summer 2009  
US Department of Energy  
GCEP-SURE program  
UALR – Department of Chemistry

# Outline

- Background : Atmospheric aerosols
- TXRF Spectroscopy
- Measuring black carbon (BC) concentration in the atmosphere and its dependence on weather conditions.
- Future work
- Acknowledgements



# Atmospheric aerosols physical and chemical properties

- Aerosols are suspension in a gas, usually air, of microscopic liquid or solid particles, such as smoke, dust, fog, or smog, that tend to remain dispersed rather than to settle.
- They help in the formation of cloud droplets acting as cloud condensation nuclei.
- More aerosol in the atmosphere equals smaller water droplets in clouds (average size 20 microns) which reflects more solar radiation back to space.



# Atmospheric aerosols physical and chemical properties

- Scattering is primarily dependent on particle size.
  - ⊙ Ex. Sulfate aerosols which cooling effect in the N Hemisphere is comparable in magnitude to the atmospheric warming produced from increases in  $\text{CO}_2$ .
- Absorption is largely dependent upon chemical composition.
  - ⊙ Ex. BC aerosols produced from incomplete combustion can reduce the amount of sunlight reaching the Earth's surface by as much as 10%.
  - ⊙ However, they can also increase the solar energy absorbed in the atmosphere by as much as 50% (dominant absorber) leading to significant local warming of the boundary layer.




# Atmospheric aerosols physical and chemical properties

- Other light absorbing carbonaceous species that are present in atmospheric aerosols are:
  - Aromatic hydrocarbons – produced along with BC during combustion
  - Nitrated derivatives – formed by atmospheric oxidation in the presence of  $\text{NO}_2$ .
  - Amino acids – biological particles
  - Polycarboxylic acids known as “humic-like” substances (HULIS) – produced by biomass burning, atmospheric oxidation of BC or atmospheric oxidation of biogenic hydrocarbons



Bio-hydro-atmosphere interactions of Energy, Aerosols, Carbon, H<sub>2</sub>O, Organics and Nitrogen



# Outline

- Background : Atmospheric aerosols
- TXRF Spectroscopy
- Measuring black carbon (BC) concentration in the atmosphere and its dependence on weather conditions.
- Future work
- Acknowledgements



# Total X-Ray Fluorescence Spectroscopy

- Bruker's TXRF S2 PICOFOX is a portable bench-top instrument that is useful for quantitative and semi-quantitative multi-element microanalysis of liquids, solids or powders.

- The principle used for such analysis is called total reflection X-ray fluorescence spectroscopy (or TXRF spectroscopy).


- This instrument can be used to trace metals in atmospheric aerosols samples collected in a filter media.






# TXRF Spectroscopy

Description of principle:



- The primary beam is generated by a X-ray tube.
- The X-ray spectrum depends on the anode material which in this case is molybdenum .
- This beam is monochromatized by Bragg-reflection on a multilayer.
- Then these X-rays hit a sample carrier under a very small angle of incidence.
- The X-rays are totally reflected by the surface of the sample carrier and fluorescence radiation is emitted by the sample substance.
- From this we get a spectrum for the sample (see next slide).

\*Sample carrier is either a quartz glass disk or an acrylic glass disk.





# TXRF PICOFOX

- All elements present in the sample are identified using the periodic table as shown in figure 1. The concentration can be obtained either in ppm or %.
- The series lines fitted for each element are shown in figure 2.
- Some advantages of this instruments are:
  - Data and/or graphic can be easily exported to an excel sheet for further analysis.
  - Sample amounts needed for study are in the nanogram to microgram range.
  - Low operating cost.




Figure 1




Figure 2




# Outline

- Background : Atmospheric aerosols
- TXRF Spectroscopy
- Measuring black carbon (BC) concentration in the atmosphere and its dependence on weather conditions.
- Future work
- Acknowledgements




# Aethalometer

- The aethalometer is a tape sampler that takes in air at a controlled flow rate, and then passes the air through a special glass fiber filter tape.
- Black carbon is measured using optical absorption at a near-IR source wavelength (880 nm).
- An additional second wavelength source is available in the near-UV (370nm)
- These two wavelengths can be used to compare emissions due to low combustion efficiency versus engine deterioration and for providing further information about different potential sources of BC.




# Aethalometer : Ex. UALR Near UV vs Near IR BC concentration measurements

24-Jul-09 – Non rainy day




5-Aug-09 – Rainy day





# Weather data collection



- Also a rain gauge is used to collect water samples during rain events.



- WXT 510 is an instrument designed to measure various weather parameters like wind speed and direction, RH, temperature, pressure, rain intensity and hail intensity.

# Conclusion

- TXRF Spectroscopy
  - ④ This type of spectroscopy has proven to be very useful to trace metals in bulk samples like the atmospheric aerosols collected in filter media.
  - ④ Also it provides the concentration in the sample of each of those metals which is very useful for characterization of compounds in the sample.
- Black carbon sources
  - ④ Studying BC concentration variation in a day using the near UV (370 nm) and N-IR (880 nm) wavelengths seems to be a good way of characterizing BC concentration diurnal cycle and see the trends during rain and non-rain events.



# Future work

- TXRF Spectroscopy
  - Run all samples available for further analysis using trace metals.
- Black carbon vs DOC concentrations
  - Study black carbon concentrations and dissolved organic carbon concentration in rain water during rain events at UALR.
  - Also correlating this amounts to weather variables such as temperature, boundary layer mixing, relative humidity and wind.



# Acknowledgments

This work was supported by the Global Change Education Program (GCEP-SURE), U.S. Department of Energy, Atmospheric Science Program.

Thanks to my mentors N.A. Marley & J.S. Gaffney that accepted me as a SURE fellow for the Summer 2009 and my co-workers for a great research experience..

Also to all UALR- Chemistry Department for their hospitality. Milton Constantin & Ms. Rose Etta Cox at Oak Ridge.

