

Kinetics of Submicron Oleic Acid Aerosols with Ozone; a Novel Aerosol Mass Spectrometric Technique

J. W. Morris and P. Davidovits

Chemistry Department, Merkert Chemistry Center, Boston College, Chestnut Hill, MA 02467

J.T. Jayne, J. L. Jimenez[°], Q. Shi, C.E. Kolb, and D.R. Worsnop

Center for Aerosol and Cloud Chemistry, Aerodyne Research Inc., 45 Manning Road, Billerica, MA 01821

W.S. Barney* and G. Cass[•]

Department of Environmental Engineering, California Institute of Technology, Pasadena, CA 91125

Abstract. The reaction kinetics of submicron oleic (9-octadecanoic (Z)-) acid aerosols with ozone was studied using a novel aerosol mass spectrometric technique. In the apparatus a flow of size-selected aerosols is introduced into a flow reactor where the particles are exposed to a known density of ozone for a controlled period of time. The aerosol flow is then directed into an aerosol mass spectrometer for particle size and composition analyses. Data from these studies were used to: (a) quantitatively model the size-dependent kinetics process, (b) determine the aerosol size change due to uptake of ozone, (c) assess reaction stoichiometry, and (d) obtain qualitative information about the volatility of the reaction products. The reactive uptake probability for ozone on oleic acid particles obtained from modeling is $1.6(-0.2) \times 10^{-3}$ with an upper limit for the reacto-diffusive length of ~ 10 nm. Atmospheric implications of the results are discussed.


[°] Present address: Dept. of Environmental Engineering, California Institute of Technology, Pasadena, CA 91125


* Present address: Arthur D. Little, Inc., Acorn Park, Cambridge MA 02140


[•] Deceased.


Laminar Flow Reactor


Aerosol Mass Spectrometer














Laboratory Studies of Sulfuric Acid Aerosol Formation, Growth, and Coagulation

J. W. Morris and P. Davidovits

Chemistry Department, Merkert Chemistry Center, Boston College, Chestnut Hill, MA 02467

John T. Jayne, Charles E. Kolb, and Douglas R. Worsnop


Center for Aerosol and Cloud Chemistry, Aerodyne Research Inc., 45 Manning Road, Billerica, MA 01821

Darren D. O'Brien and Gregory J. McRae

Department of Chemical Engineering, Massachusetts Institute of Technology, Cambridge, MA 02139

Abstract. Laboratory measurements involving nucleation, condensation, and coagulation of sulfuric acid are described. Studies were carried out using a laminar flowtube at atmospheric pressure coupled to an aerosol mass spectrometer (AMS). A sulfuric acid reservoir heated to a known temperature was connected to the flowtube via heated tubing. Dry filtered nitrogen carrier gas was pre-heated and passed over the reservoir entraining an equilibrium density of the $\text{H}_2\text{SO}_4(\text{g})$ to the flowtube. The entrance to the flowtube was heated to a temperature above that of the reservoir, establishing a well-characterized decreasing temperature profile along the flowtube axis. Supersaturation (and consequent nucleation and growth) conditions were reached at a position within the flowtube determined by the temperature profile and the H_2SO_4 vapor pressure. The number and mass distributions of the particles formed in this way were measured on-line by the AMS at the flowtube exit. With the entrance of the flowtube at 150 C, particles were first observed at a H_2SO_4 vapor pressure of 10^{-4} torr forming particles at a number density of 10^4 cm^{-3} . As the pressure increased to 3×10^{-2} torr, the particle number density increased to 10^7 cm^{-3} . However, the median diameter increased only from 40 to 50 nm. At pressures higher than 4×10^{-3} torr, the size distributions were observed to exhibit broadening and a new mode appeared near 60nm, consistent with coagulation. A third mode appeared at H_2SO_4 pressures near 3×10^{-2} torr indicating further coagulation. This work is a first time *on-line* laboratory observation of sulfuric acid microphysics spanning nucleation and coagulation processes. The features of the apparatus enabling absolute detection of mass and number distributions will be discussed.


Experimental Apparatus



Aerosol
Mass
Spectrometer



~150 C




*Particle formation,
coagulation,
and growth*

*Particle
Detection*

$T_{H_2SO_4} < 150\text{ C}$

H₂SO₄ Entrainment



20x10³

15

10

5

0

10

88 to 92 C

2

3

4

5

6

7

8

9

100

2

3

4

5

6

7

8



